

Cortical cancellous chips living donor (Frozen)

		Dimensions	
		Volume	Range
Description		g	mm
T0629-05	Cortical cancellous chips	approx. 5 g	<=4 mm
T0629-30	Cortical cancellous chips	approx. 30 g	<=4 mm
T0630-15	Cortical cancellous chips	approx. 15 g	<= 8mm
T0630-30	Cortical cancellous chips	approx. 30 g	<= 8mm
T0628-30	Cortical cancellous chips	approx. 30 g	<= 12mm

Processing specifications

Femoral head from living donor; aseptically cleaned and milled (Spierings milling drum); bone chips $5 \times 8 \times 10^{-5}$ x washed with sterile saline. Composed of approximately 85% cancellous bone (+/-5%) with 95% of the pieces ranging in size per product code.

Measurements

Weight weight measured in gram
Particle size measured in millimeters

FG -T0628/T0629/T0630

© 2018