



DBM - Putty

Description

621005	DBM – Putty	0.5 cc
621010	DBM – Putty	1 cc
621020	DBM – Putty	2 cc
621025/BT1013	DBM – Putty	2.5 cc
621050	DBM – Putty	5 cc

Processing specifications

DBM-Putty contains demineralized cortical bone powder and a carrier. Putty has a consistency similar to a gel or mastic and is ready for use. Demineralized powder of the putty was found to be osteoinductive in an animal model.

Measurements

V volume measured in cc