



Femoral shaft

Dimensions

Description

	<i>L</i>	<i>(W)</i>	<i>(A)</i>	<i>(D')</i>	<i>(L)</i>
	Cm	mm	mm	mm	cm
550430 Femoral shaft 10.0 cm	9.9-10.1	10.0-50.0	10.0-50.0	5.0-30.0	10.0
550450 Femoral shaft 20.0 cm	19.9-20.0	10.0-50.0	10.0-50.0	5.0-30.0	20.0
550455 Femoral Shaft > 20.0 cm	20.1-50.0	10.0-50.0	10.0-50.0	5.0-30.0	20.1-50.0

Processing Specifications

Cortical wall thickness must be 3mm or greater at all points

Measurements

- L overall length
- W width of the shaft measured anterior to posterior
- A width of the shaft measured medial to lateral
- D' Diameter of the medullary canal measured medial to lateral

FG-205
©2018

