



Bone Tendon Bone

Dimensions

Description

		(L)	(W)	(X)	(Z)	(Q)	(y)	(T)
		cm	mm	mm	mm	mm	mm	mm
530005	With quadriceps, hemi bone block and small tendon	8.0-14.0	10-13	30-75	07-13	30-80	25-30	10-13
530010	With quadriceps, hemi bone block	8.0-14.0	13-30	30-72	07-25	30-80	25-30	16-30
530030	Whole with quadriceps 50-80 mm	13.0-20.0	13-40	30-80	07-30	50-80	80-85	16-30
530036	Whole with quadriceps >80 mm	13.0-20.0	13-40	30-80	07-30	>80	80-85	16-30
530031	Whole without quadriceps	13.0-20.0	13-40	30-80	07-30	N/A	80-85	16-30

Processing Specifications

Tendon is attached to patellar block at the proximal end and tibial block at the distal end.


Measurements

All measurements except length are recorded in millimeters. Length is recorded in tenths of a centimeter

L=length of the tendon measured along the patellar and tibial blocks

T=thickness of the tibial block measured at the minimum point anterior to posterior

W=minimum medial to lateral width of the tendon



X=distance from the proximal tibia and distal patellar block

Y=length of the tibial block

Z=minimum width of the tibial block medial to lateral

Q=length of the quadriceps tendon from attachment on the superior patella to the cut end

Illustration:

