



Tibia proximal with BTB

Dimensions

Description

		(L)	(W)	(D')	(A')	Q	(W')	(X)	(W1)	(T1)	(T2)
		cm	mm	mm	mm	mm	mm	mm	mm	mm	mm
550060	Tibia proximal Left with BTB	20-99	10-50	1-30	10-50	50-200	30-100	30-80	13.4-40	1-200	1-200
550070	Tibia proximal Right with BTB	20-99	10-50	1-30	10-50	50-200	30-100	30-80	13.4-40	1-200	1-200


Processing specifications

Quadriceps tendon must be intact and no less than 50 mm in length.

Measurements

L	Maximum length from the top of the tibial plateau to the cut end
W	Width of the shaft at cut end measured anterior to posterior
D'	Diameter of the medullary canal measured medial to lateral.
A'	Width of the shaft at cut end measured medial to lateral
Q	Length of the quadriceps tendon
W'	Maximum width of the plateau measured medial to lateral
X	Minimum length of the patellar tendon; distance from insertion site on tibia to insertion site on patella
W1	Minimum medial to lateral width of the patellar tendon

- T1 Medial condyle measured anterior to posterior
- T2 Lateral condyle measured anterior to posterior



PROXIMAL